



BIONEK 

Hospital Furniture

One Function Manual Bed

Model: BX-A120B

1. Bed size: 2040*925*500mm
2. Mattress platform: 1950*810mm
3. Backrest adjustment: 0-70 ° ± 5 °
4. Head and foot board: PE material shaped by blow molding- retractable
5. Bed surface: PUNCHING PRESS BED BOARD with powder coated
6. Six-column aluminum alloy folding side rails
7. 125mm casters with individual brakes, portable, retractable
8. Accessories:
 - I.V pole sockets: 4pcs
 - Drainage bag hooks: 4pcs
 - IV pole (S.S IV pole, 4 hooks, height can be adjusted)
 - Mattress holder: 1pc



Two Function Manual Bed

Model: BX-A131C

1. Bed size: 2040*925*500mm
2. Mattress platform: 1950*810mm
3. Backrest adjustment: 0-70 ° ± 5 °
4. Footrest adjustment: 0-40 ° ± 5 °
5. Head and foot board: PE material shaped by blow molding- retractable
6. Bed surface: PUNCHING PRESS BED BOARD with powder coated
7. Six-column aluminum alloy folding side rails
8. 125mm casters with individual brakes, portable, retractable
9. Accessories:
 - I.V pole sockets: 4pcs
 - Drainage bag hooks: 4pcs
 - Mattress holder: 1pc
 - Protective angle: 4 pcs



3 Cranks Manual Bed

Model: BX-A131A

1. Bed size: 2150*980*500-750mm
2. Backrest adjustment: 0-70 ° ± 5 °
3. Footrest adjustment: 0-40 ° ± 5 °
4. Height adjustment: 500-750mm
5. Head and foot board: plastic shaped by blow molding- retractable
6. Bed surface: PUNCHING PRESS BED BOARD with powder coated
7. Six-column aluminum alloy folding side rails
8. 125mm casters with individual brakes, portable, retractable
9. Accessories:
 - I.V pole sockets: 4pcs
 - Drainage bag hooks: 2pcs
 - Mattress holder: 1pc



BIONEX
Rua Mario Moreira, 2675-660
Odivelas, Lisboa, Portugal
Mob: +351920096274
email: info@bionex.pt
web: www.bionex.pt