



BX-B100L

Autoclave

**CENTRAL STERILE
SUPPLY DEPARTMENT**





Features:

- Full SUS 304 stainless steel
- Accurate double scale indication pressure gauge
- Silicone seal ensuring no steam leakage
- Spin buttons to set sterilization temperature and time
- Blue, orange and red light showing work status
- Automatic release when overpressure
- Automatic power off when water shortage
- Drying function can be added

Technical Specifications:

Model	BX-B35L	BX-B50L	BX-B75L	BX-B100L	BX-B120L	BX-B150L
Volume	35L	50L	75L	100L	120L	150L
Max Working Temperature	134°C	134°C	134°C	134°C	134°C	134°C
Max Working Pressure	0.22Mpa	0.22Mpa	0.22Mpa	0.22Mpa	0.22Mpa	0.22Mpa
Time Range	0-80min	0-80min	0-80min	0-80min	0-80min	0-80min
Temperature Range	115°C-129°C	115°C-129°C	115°C-129°C	115°C-129°C	115°C-129°C	115°C-129°C
Sterilization Chamber	Φ320*485 mm	Φ320*545 mm	Φ388*480 mm	Φ388*630 mm	Φ460*590 mm	Φ460*700 mm
Power	AC 220V/50Hz/ 2.5KW	AC 220V/50Hz/ 2.5KW	AC 220V/50Hz/ 3.5KW	AC 220V/50Hz/ 4.0KW	AC 220V/50Hz/ 4.0KW	AC 220V/50Hz/ 4.0KW



BIONEX
Rua Mario Moraira, 2675-660
Odivelas, Lisboa, Portugal
Mob: +351920096274
email: info@bionex.pt
web: www.bionex.pt